



1. Do not withhold defibrillation based on weight and size criteria. Recommended settings are 2J/kg for first shock, then 4J/kg for subsequent shocks. If unable to adjust Energy settings (AED) use 200J Biphasic or Monophasic.
2. Once ETT placed, give continuous chest compressions without pauses for breaths. Give 8 - 10 "breaths" per minute.
3. Epinephrine **IV/IO**: .1cc/kg of 1:10,000 (.01mg/kg) to a maximum per dose of 10cc (1mg). May use same dose **IO** for children less than 6 years old. **ETT**: .1cc/kg of 1:1,000 (.1mg/kg) to a maximum per dose of 10mg.
4. Give 2 mins of CPR between defibrillations.